Lucitone Digital Fit™ Denture Disc Prothesenbasis Ronde Disque de dentier Disco per protesi Disco de dentadura

Instructions for Use
Gebrauchsanweisung
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso

 Manufacturer:
Dentsply Sirona
1301 Smile Way
York, PA 17404 U.S.A.
Tel. 800-243-1942
dentsplysirona.com

 **EC REP**
DeguDent GmbH
Rodenbacher Chaussee 4
63457 Hanau-Wolfgang
Germany
Tel. +49 (0) 6181-5950

See website for
symbols glossary
 
2797
DP-0000371 Rev. 9 (12/2021)

THE DENTAL
SOLUTIONS
COMPANY™

Indications

Lucitone Digital Fit Denture Discs are intended for the fabrication of PMMA based acrylic dental prostheses including but not limited to partial and full dentures.

en

Contraindications

Lucitone Digital Fit Denture Discs are contraindicated for patients who have a known hypersensitivity or severe allergic reactions to any of the components.

Warnings

Lucitone Digital Fit Denture Disc contains residual methyl methacrylate monomer which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice. Consult applicable SDS (Safety Data Sheets) for safe handling and usage. Use protective equipment to guard against dust when grinding.

Adverse Reactions

Particulates will be generated when milling or grinding acrylate resins. Safety glasses and / or gloves are recommended when milling the discs to minimize the potential for eye, skin and respiratory irritation. Adverse effects from use of this medical product are highly unlikely with proper handling and processing. However, reactions of the immune system (such as allergies) to substances contained in the material or localized paraesthesia (such as taste disturbances or irritation of the oral mucosa) cannot be completely ruled out.

Precautions

The following safety and warning information describes how to use the Lucitone Digital Fit Denture Disc in a safe manner. Refer to the special notes throughout the instructions for best practices for using the Lucitone Digital Fit Denture Disc.

- Keep product dust away from eyes.
- Do not inhale dust particles during grinding.
- Use local vacuum suction and suitable mouth/face protection during manual machining at the workplace.
- If ingested contact regional Poison Control Center.
- Comply with pertinent SDS (safety data sheets)
- Do not mill the discs with non-compatible CAD/CAM systems.
- For a list of verified CAD-CAM systems, contact Dentsply Sirona
- Follow appropriate instructions for processing of denture teeth and bonding agents.
- Users with sensitivities, dermal issues, cuts, or abrasions should wear protective gloves.
- When grinding resin materials, proper ventilation, masks, and vacuum systems should be used.

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

Workflow Solution

The Lucitone Digital Fit Denture Discs are part of Dentsply Sirona's Digital Denture Workflow Solution, which includes:

- inEos X5 Scanning
- inLab® MC X5 Milling
- Lucitone® HIPA High Impact Pour Acrylic
- Dentsply Sirona Denture Teeth (IPN 3D™, Portrait® IPN® or Genios®)

Step by Step Instructions

Lucitone Digital Fit Denture Discs are recommended for the production of full or partial denture upper and lower arches for bonded digital dentures. The discs are available in 98.5mm diameter. This dimension will fit the most popular 5-axis milling machines.

- 1. General Case Preparation:**
 - Receive and process records from dentists, following the dental laboratory procedures.
- 2. Scan and Design:**
 - Scan records following the equipment manufacturer's scanning instructions.
 - Digitally design denture following design criteria in the CAD-CAM Software.
- 3. Mill and Bond Teeth:**
 - Select the appropriate Lucitone Digital Fit Denture Disc shade.
 - Mount the disc in the milling machine's jig following the equipment manufacturer's instructions.
 - After milling, remove the disc from the machine.
 - Following the Dentsply Sirona Denture Tooth Bonding Guide for Digital Dentures, bond Dentsply Sirona denture teeth using Lucitone HIPA High Impact Pour Acrylic.
 - Remount the disc in the milling machine for a two-step mill process.
 - After final milling, remove the disc from the machine and remove the sprues.
- 4. Finish and Polish:**
 - Follow your laboratory's recommended technique for finishing, festooning, and polishing dentures.
- 5. Repair and Reline:**
 - Repairs and relines of Lucitone Digital Fit dentures can be carried out with conventional auto-polymerizing ('cold cure') acrylics. Dentsply Sirona recommends Lucitone HIPA High Impact Pour Acrylic. Follow the manufacturer's instruction for the material used.

Important: Dentsply Sirona is committed to providing open system solution support for dental laboratory production of Digital Dentures. Fundamental operability has been tested and pre-qualified with the list of "open" milling machine systems. Contact your Dentsply Sirona representative for a comprehensive list of Digital Denture machine prerequisites and supported systems. For assistance regarding specific machine operability relating to digital dentures, preventive maintenance, service, and particularly related to firmware and CAM software, please contact your machine reseller or manufacturer.

DENTSPLY SIRONA HAS TESTED AND PRE-QUALIFIED A LIST OF OPEN MILLING MACHINE SYSTEMS. DENTSPLY SIRONA ONLY RECOMMENDS USE OF THESE SYSTEMS. DENTSPLY SIRONA SHALL NOT HAVE ANY LIABILITY TO THE DENTAL LABORATORY OR ANY OTHER PERSON OR ENTITY FOR ANY INDIRECT, INCIDENTAL, SPECIAL, OR CONSEQUENTIAL DAMAGES WHATSOEVER, INCLUDING, BUT NOT LIMITED TO, LOSS OF REVENUE OR PROFIT, LOST OR DAMAGED DATA OR OTHER COMMERCIAL OR ECONOMIC LOSS, EVEN IF DENTSPLY SIRONA HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, OR THEY ARE FORESEEABLE. DENTSPLY SIRONA AND IT'S AFFILIATES ARE ALSO NOT RESPONSIBLE FOR AND SHALL NOT HAVE ANY OBLIGATION OR LIABILITY FOR CLAIMS BY A THIRD PARTY.

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