OptraSculpt® & OptraSculpt® Pad

[en] Instructions for Use
Composite modelling instruments

Rx ONLY

CE 0123

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

Date information prepared:
2024-01-25 / WW / Rev. 2
web

ivoclar

English

1 Intended use

Intended purpose

Shaping and contouring of uncured sculptable composites

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth
- Adult patients with dental implants
- Edentulous adult patients

Intended users / Special training

- Dentists (clinical procedure)
- Dental laboratory technicians (shaping and smoothing in the dental laboratory)

Application

For dental use only.

Description

OptraSculpt® and OptraSculpt® Pad are modelling instruments with corresponding attachments designed for adapting, shaping, contouring and smoothing uncured, sculptable direct composite resins in direct restorative procedures as well as for the modelling and smoothing of laboratory composites in paste form.

Indications

None

Areas of application

- OptraSculpt is suitable for adapting, shaping and contouring uncured composites.
- OptraSculpt Pad is suitable for smoothing and adapting uncured composites.

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If the stipulated working technique cannot be applied.
-  The modelling attachments are intended for single use only. They must not be cleaned, disinfected, sterilized and reused.

Processing restrictions

- Cleaning and disinfection agents can negatively affect the modelling attachments, autoclaving will destroy the attachments.
- The hand instrument is autoclavable and suitable for multiple use. It is designed to withstand at least 200 autoclaving cycles. Before use, visually inspect the product for damage.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

OptraSculpt:

- OptraSculpt attachments: Polypropylene and styrene ethylene butylene styrene
- OptraSculpt hand instrument: Stainless steel and polyether ether ketone

OptraSculpt Pad:

- OptraSculpt Pad attachments: Polyethylene and ethylene vinyl acetate
- OptraSculpt Pad hand instrument: Stainless steel and polyether ether ketone

2 Application


Select the appropriate OptraSculpt or OptraSculpt Pad attachment, depending on the expected size of the filling and the area in which it will be employed. Remove the OptraSculpt Pad attachment from the packaging indirectly (e.g. using sterile tweezers). Then insert the attachment into the ring-shaped end of the hand instrument. The colour-coded OptraSculpt attachments are supplied correctly aligned in refill boxes, so that they can be removed and taken up directly with the ring-shaped end of the hand instrument. To attach the disc-shaped OptraSculpt Pad attachments, place them flat on a firm, contamination-free working surface and insert them into the ring-shaped end of the hand instrument.

Place the disc-shaped OptraSculpt Pad attachment flat and upright on a working surface.

Open the dispenser box by exerting slight pressure on the marked surface.

Insert the ring over the body until it meets the base.

Select the appropriate attachment based on the colour coding described below:
Anthracite: ball; Turquoise: chisel; Purple: point

OptraSculpt Pad attachment is correctly inserted in the holding ring.

Slide the ring over the shaft of the attachment until the ring sits flat on the attachment. Then remove the attachment, which is now firmly attached to the instrument.


The OptraSculpt and OptraSculpt Pad modelling attachments are intended for single use only. After having applied the uncured, sculptable composite restorative, adapt and shape it with a dabbing motion, applying gentle pressure. After intraoral use and prior to disinfecting, cleaning and sterilizing the hand instrument, OptraSculpt and OptraSculpt Pad attachments must be removed and disposed of appropriately.

Using the reference scales on the OptraSculpt Pad instrument handle:

The OptraSculpt Pad hand instrument features two reference scales.

Scale 1 provides plain reference markings to compare the tooth-to-tooth width proportions. The horizontal line serves as a reference to evaluate the course of the incisal edge line.


Scale 1 can be aligned with the incisal edge to e.g. compare the contralateral tooth width proportions or evaluate the course of the incisal edge line against the horizontal line on the scale.

Scale 2 contains the markings for the ideal average anterior tooth widths (dashed vertical lines) and the typical average position of the tooth axes in the upper jaw (oblique dotted lines). They provide esthetic orientation in the realization of extensive anterior reconstructions. The values are relative values that are intended to serve as a guide in the harmonious, anatomical design of anterior restorations. The scale is positioned correctly if the triangle in the middle of the scale points towards the contact point between the central incisors in the upper jaw and the numbers at both sides of the scale are consistent with the quadrants under examination.

Scale 2 is unsuitable for use in the mandibular jaw!

Scale 2 is utilized to compare the clinical situation against the ideal average tooth width proportion and angular alignments in the maxillary anterior dentition. Consequently, esthetic enhancements can be easily planned and areas of improvement identified.

Important:

This scale indicates relative values only. The marked widths and angular alignments should not simply be taken at face value.

Example: The example depicted above shows comparatively ideal tooth width proportions, even if the overall tooth width may appear slightly larger than average. The angles of the lateral incisors and canines, however, are too straight.

Application in the dental laboratory

OptraSculpt and OptraSculpt Pad are also suitable for the shaping and contouring of laboratory composites in paste form, e.g. SR Nexco® Paste. Due to the attachment's non-stick properties, large veneers can be conveniently contoured and given an ideal shape, particularly with the small pad. OptraSculpt and OptraSculpt Pad are a useful addition to the range of instruments in the dental laboratory. They are especially useful for the processing of laboratory composites in paste form.

3 Cleaning, disinfection and sterilization

General notes

- Please observe the country-specific regulations and guidelines regarding the hygiene standards and reprocessing of medical devices in dental practices.
- Clean and disinfect the hand instruments before each use.
- The hand instrument must be sterilized immediately before each use.
- Immediately after each clinical application, the hand instrument must be cleaned and disinfected.
- Do not place the instruments in NaCl solution (to avoid the risk of pitting and/or stress corrosion cracking).
- Protective gloves and goggles complying with the requirements of the regulation on personal protective equipment PPE (EU) 2016/425 must be worn when handling any used and contaminated instruments.

- A washer-disinfector (according to EN ISO 15883) whose effectiveness has been verified must be used to ensure reliable machine cleaning and disinfection. It is the responsibility of the operator to define, document and implement the validation, continued compliance with the performance requirements and regular routine control of the reprocessing performance of the washer-disinfector.

Cleaning preparation

- Carefully remove large unhardened residues from the instruments using a lint-free cellulose tissue. For preliminary cleaning, place the instruments in a disinfection bath (with a lid and sieve insert) filled with aldehyde-free alkaline disinfectant (e.g. immersion in 2% ID 212 forte, Dürr Dental, in tap water / drinking water quality, for 15 min). Make sure that the instruments are completely covered in disinfectant and that the disinfectant is free of bubbles.
- Only use a soft brush for the manual removal of residues and brush each instrument for at least 20 seconds until no residues are visible. Never use a metal brush or steel wool.
- Rinse with distilled water for 1 min.
- After preliminary cleaning, carefully rinse the instruments under running tap water (drinking water quality) for at least 10 seconds per instrument.

Automated cleaning and disinfection

- Thermal disinfectant: Manufacturer's specification according to DIN EN ISO 15883. Cleaning program as indicated by the manufacturer's operating instructions.
- Place the instruments in a sieve suitable for small parts, or position them on the loading rack in a way that ensures cleaning and disinfection of all surfaces of the instruments.
- Clean and disinfect the instruments in a washer-disinfector using a thermal program (A_0 value > 3000 is achieved at 90 °C/194 °F for 5 minutes) according to the directions of the device manufacturer.
- Cold water rinse for 1 min
- As cleaning agent, use e.g. neodischer FA (0.2% v/v) (5 min, 55 °C / 131 °F)
- Neutralization (if required) with e.g. neodischer Z (0.1% v/v)
- Rinse with distilled water for 1 min.
- Subsequently, dry the instruments according to the directions of the device manufacturer.
- Only use deionized water with a low germ level (max. 10 cfu/ml) and a low endotoxin level (max. 0.25 Eu/ml) for final rinsing.

Manual cleaning and disinfection

Cleaning

- Before applying the manual cleaning procedure described below, any gross contamination (tissue or bone residues, etc.) must be removed by means of suitable blood fixation methods.
- The instruments should be immersed in a suitable cleaning solution according to the directions of the cleaning agent manufacturer, e.g. 2.0% (v/v) ID 212 forte in tap water (minimum quality: drinking water; temperature: 20 ± 2 °C/68 ± 3.6 °F) for 15 minutes. Make sure that the immersed surfaces are completely covered in cleaning agent and that the cleaning agent is free of bubbles.
- Then clean the instruments manually by thoroughly brushing with a soft, non-metal brush until no residues are visible. Subsequently rinse the instruments under running tap water (minimum quality: drinking water; temperature: 20 ± 2 °C/68 ± 3.6 °F) for 10 seconds each.

Disinfection

- Immerse the instruments in a suitable disinfectant for 15 min at 20 °C ± 2/68 ± 3.6 °F (e.g. Dürr ID 212 forte; 2% v/v in drinking water). Make sure that the immersed instruments are completely covered in disinfectant and that the disinfectant is free of bubbles.
- Remove any remaining disinfectant according to the directions of the manufacturer of the disinfectant (e.g. rinsing with water).
- Only use deionized water with a low germ level (max. 10 cfu/ml) and endotoxin level (max. 0.25 Eu/ml) for final rinsing.

Drying

Compressed air or clean, lint-free cellulose tissue

Sterilization

Instruments must be sterilized immediately before use. Ivoclar Vivadent AG recommends one of the following sterilization procedures:

- For countries outside the United States:
To pack the items for sterilization, only use a sterile barrier system that complies with ISO 11607-1 (e.g. steriCLIN) made of paper/film or an aluminium sterilization container, which is designated by the manufacturer for steam sterilization. The sterile barrier system being used must be large enough. The filled sterile barrier system must not be stretched.
- For the United States:
Insert the products in a perforated basket with lid and wrap in two layers of 1-ply polypropylene wrap using sequential envelope folding technique prior to sterilization.
Note: Users in the United States must make sure that the sterilizer and any sterilization accessories (e.g. sterilization wraps, pouches, basket, biological or chemical indicators) are approved by the FDA for the intended sterilization.

Steam-sterilize using a fractionated pre-vacuum process in a steam sterilizer (e.g. Systec HX-320 autoclave) according to ISO 17665 under the following conditions:

	Method	Conditions	Drying time
1	Steam sterilization (autoclave) Fractionated vacuum	134 °C (273 °F) for 4 min	Local practice
2	Steam sterilization (autoclave)* Fractionated vacuum	132 °C (270 °F) for 3 min	10 min
3	Steam sterilization (autoclave)** Fractionated vacuum	134 °C (273 °F) for 3 min	Local practice
4	Steam sterilization (autoclave)*** Fractionated vacuum	134 °C (273 °F) for 18 min	Local practice

[*] recommended for the USA

[**] recommended for the UK

[***] recommended for Switzerland and France

Note: No sterilization in chemiclave or hot air disinfectors!

Storage

Sterilized products packed in a sterile barrier system (e.g. sterilization pouch) are intended for immediate use and must not be stored for longer than 48 hours.

Maintenance, examination and inspection

After cleaning and disinfecting the instruments, visually inspect them to see if they are clean, undamaged and function properly. Repeat the reprocessing procedure if there is evidence of any macroscopically visible residual contamination. If any external defects are detected, do not use the instrument any more.

4 Safety information

- Before autoclaving the hand instrument, remove the modeling attachments.
- Forceful removal of the attachments or mechanical manipulation of the holding ring may lead to instrument breakage and should therefore be avoided.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available on the website www.ivoclar.com
- Explanation of the symbols: www.ivoclar.com/eIFU

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- It is recommended to use a rubber dam to avoid accidental ingestion or aspiration of small parts.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:
Ingestion or aspiration of material (intraoral use)

5 Shelf life and storage

- Store the instruments in a dry place to avoid condensation build-up.
- Store OptraSculpt and OptraSculpt Pad attachments in the original packaging at room temperature, protected from dust, moisture, compression and contamination.
- Storage temperature 2-28 °C (36-82 °F)

Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

6 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

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