# X-Smart<sup>®</sup>

## Apex Locator

EN	USER MANUAL .....	5
FR	MANUEL D'UTILISATION .....	43
IT	MANUALE D'USO .....	83
DE	BENUTZERHANDBUCH .....	123
ES	MANUAL DE USUARIO .....	163
PT	MANUAL DO UTILIZADOR .....	203
NL	GEBRUIKERSHANDLEIDING .....	243
SV	BRUKSANVISNING .....	283
DA	BRUGERMANUAL .....	323
NO	BRUKERHÅNDBOK .....	363
FI	KÄYTTÖOHJE .....	403
PL	INSTRUKCJA OBSŁUGI .....	443
CS	UZIVATELSKA PRIRUČKA .....	484
HU	FELHASZNALÓI KEZIKÖNYV .....	525

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PAGE INTENTIONALLY LEFT BLANK

## Table of contents

<b>Introduction</b> .....	<b>7</b>
<b>1 Indication For Use</b> .....	<b>8</b>
<b>2 Contraindications</b> .....	<b>8</b>
<b>3 Warnings</b> .....	<b>8</b>
<b>4 Precautions</b> .....	<b>9</b>
<b>5 Adverse Reactions</b> .....	<b>10</b>
<b>6 Step-by-Step Instructions</b> .....	<b>10</b>
6.1 Content.....	10
6.2 X-Smart® Apex Locator Overview .....	11
6.3 Connecting the AC Plug Adapter .....	12
6.4 Recharging the Battery .....	12
6.5 Cable Connection Test .....	13
6.6 Apex Localization .....	14
6.6.1 Getting Started .....	14
6.6.2 Apex Localization .....	15
6.6.3 Over-Instrumentation .....	16
6.6.4 Completion of the measurements .....	16
6.7 Sound Adjustment.....	16
6.8 Device Settings .....	17
6.8.1 Dr's Choice.....	18
6.8.2 Language selection.....	19
6.8.3 Check Mode .....	19
6.8.4 Demo Mode.....	20
6.8.5 Sleep Mode .....	20
6.8.6 Automatic Shutdown .....	21
6.8.7 ECO Mode .....	21

<b>7</b>	<b>Maintenance, Cleaning, Disinfection and Sterilization .....</b>	<b>22</b>
7.1	General Recommendations .....	22
7.2	Disinfection and Sterilization Procedure .....	23
<b>8</b>	<b>Troubleshooting.....</b>	<b>26</b>
<b>9</b>	<b>Warranty.....</b>	<b>28</b>
<b>10</b>	<b>Disclaimer.....</b>	<b>29</b>
<b>11</b>	<b>Certification .....</b>	<b>29</b>
<b>12</b>	<b>European Authorized Representative.....</b>	<b>30</b>
<b>13</b>	<b>Disposal of the Product.....</b>	<b>30</b>
<b>14</b>	<b>Reporting of Incident to Manufacturer .....</b>	<b>30</b>
<b>15</b>	<b>Technical Characteristics .....</b>	<b>30</b>
<b>16</b>	<b>Identification of Symbols .....</b>	<b>32</b>
	<b>Appendix.....</b>	<b>34</b>

**For additional languages, visit our website:  
[dentsplysirona.com/ifu](https://dentsplysirona.com/ifu)**

Technical modifications on our product are not subject to notification.  
Photos of our devices are not contractual.

## Introduction

Congratulations on the purchase of **X-Smart® Apex Locator**.

**X-Smart® Apex Locator** is indicated for root canal treatment, which requires the precise determination of position of the endodontic file in the canal with respect to the root apex. The device is aimed at detecting the minor apical foramen based on analysis of electrical properties of different tissues inside the root canal system.

**X-Smart® Apex Locator** enables precise determination of the apical position of the endodontic file tip, providing correct results in canals with various conditions – dry and wet. The device is designed for a high level of health and safety of patients and users. The clinical benefits include better precision and confidence of the dentist and x-ray exposure reduction for the patients.

For optimal safety and performance, read this user manual carefully before use. Make sure you have understood and followed the clinical precautions - as well as the general warnings, precautions and contraindications - before proceeding to determine a working length. Keep this user manual for future reference.

# 1 Indication For Use

**X-Smart® Apex Locator** is an electronic device used for apex localization and working length determination during root canal treatment.

The intended patient population consists of patients who need to undergo root canal treatment.

# 2 Contraindications

**X-Smart® Apex Locator** is not recommended for use in patients that have a pacemaker or other implanted electrical devices.

# 3 Warnings

- **X-Smart® Apex Locator** must only be used in hospital environments, clinics or dental offices either by skilled endodontic experts or by qualified dental practitioners performing root canal treatments.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of **X-Smart® Apex Locator**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

## 4 Precautions

- Do not use **X-Smart® Apex Locator** in the vicinity of devices emitting electromagnetic noise such as x-ray viewer with fluorescent lamps, film viewers, ultrasonic devices, etc.
- During device operation protect **X-Smart® Apex Locator** from occasional spillage of liquids.
- Do not use **X-Smart® Apex Locator** in presence of flammable materials.
- **X-Smart® Apex Locator** should be used with the manufacturer's original accessories only.
- In order to prevent infectious agent transmission it is highly recommended to use a rubber dam system during the endodontic procedure.
- To ensure that short circuits do not impair the measurements, be particularly careful with patients fitted with metallic crowns, bridges or large metallic fillings (avoid any contact of the file or the lip clip with metals).
- High concentrations of sodium hypochlorite may result in a lower accuracy of the measurements. For working length determination, we recommend to use sodium hypochlorite solution at maximum 5% concentration.
- Make sure that the canal is wet enough to ensure reliability of the measurement.
- Ensure that the file does not touch other instruments.
- Avoid excessive liquids inside the tooth cavity to prevent overflow and incorrect measurements.
- Teeth with open apices may give imprecise results.
- The use of apex locators alone without a preoperative and postoperative radiograph is not a recommended practice, since apex locators may not be able to work properly in all conditions. It is mandatory to confirm radiographically the working length established using the apex locator.
- For your own safety, please use personal protection gear (gloves, mask).

## 5 Adverse Reactions


None.

## 6 Step-by-Step Instructions

### 6.1 Content

Check the content of the equipment before use:


- **X-Smart® Apex Locator**
- Charger
- Measurement cable
- Lip clip (2 pcs)
- File clip (2 pcs)
- Touch probe
- User manual
- Tip card


Lip clip  
(2 pcs.)

Touch probe

File clip  
(2 pcs.)

Measurement  
cable

#### Note

Measurement cable with attached lip clip and file clip (or touch probe) constitute Applied Parts of the device.

---

## 6.2 X-Smart® Apex Locator Overview


Fig. 1 Front view


Fig. 2 Rear view

## 6.3 Connecting the AC Plug Adapter

Select the plug adapter that matches your electric power outlet.


Fig. 3 Plug adapters for the charger

Slide the plug adapter downwards into the slots until it locks in place with a click. To remove, press the locking button (see Fig. 3) and pull out the plug adapter.

## 6.4 Recharging the Battery

**X-Smart® Apex Locator** is a battery-operated portable device which is powered by a lithium-ion rechargeable battery. Battery status during the operation is shown on the screen.

When the battery is low (lower than 25%), battery indicator will turn blinking red, meaning that the battery requires recharging. However, the **X-Smart® Apex Locator** will continue normal operation even with low battery for several treatments before the device shuts down.

To charge the battery, follow the next steps:

1. Complete the measurements and disconnect the measurement cable from the patient.
2. Unplug the measurement cable from the device.


3. Connect the charger cable to **X-Smart® Apex Locator** USB-C receptacle located on the right side of the device. See Fig. 1
4. Connect the charger to the mains.

### Warning

While charging, the charger and the device should be outside patient environment (at least 1.5m from the patient).

During battery charging, charging screen is activated with battery symbol and battery level (percentage) on the screen.

Duration of charging: Approximately 3 hours.

### Notes

- Use only the original charger.
- **X-Smart® Apex Locator** cannot be used while charging.

In case that the battery is completely flat, and the device would not turn on, contact your distributor for battery replacement by a service personnel.

## 6.5 Cable Connection Test

In case measurement indication is not detected a cable connection test needs to be conducted. A connection test feature is included in **X-Smart® Apex Locator** in order to check the cables:

1. Connect the measurement cable with attached file clip and lip clip and turn on the device.
2. Connect the metal part of the file clip to the lip clip. Make sure that the accessories are cleaned properly before the test.
3. “Cable connection test” icon  
 should appear on the screen - see Fig. 4.
4. If no icon appears, the file clip or the measurement cable should be replaced.


Fig. 4

## 6.6 Apex Localization

### 6.6.1 Getting Started

Disconnect the charger from the device if connected.

1. Before connecting the measurement cable with attached lip clip and file clip to the patient, plug the measurement cable into the device and turn on the device by pressing the “ON / OFF” button.
2. Attach the lip clip to the patient – see Fig. 5.
3. Gently insert the file into the canal and connect the file clip to the metal shaft of the file – see Fig 6. Alternatively, use the touch probe to provide electrical contact with the endodontic file – see Fig. 7.


Fig. 5


Fig. 6


Fig. 7

#### Note

To ensure optimal performance the file size should be adjusted to the canal diameter.

4. Two initial beeps indicate closed measurement circuit and beginning of length determination. The file movement in the canal is shown on the full canal image on the screen.


Fig. 8

**Note**

Absence of the two beeps audio signal indicates a faulty connection. Disconnect the measurement cable from the patient and check cable connections, clean the file clip and the lip clip, moisten the canal if necessary and start again.

No other adjustments are required before starting apex localization.

**6.6.2 Apex Localization**

Advance the endodontic file slowly into the canal. The first three blue bars correspond to the coronal/medial section. When the file reaches the fourth blue bar, very close to the apical region, an audio signal sounds. As the file progresses in the canal, subsequent green bars in the apical region turn on gradually and the interval between the audio signals becomes shorter – see Fig. 8.

**Note**

The bars indication on the **X-Smart® Apex Locator** screen does not represent a distinct length or distance in mm or other linear units. It simply indicates the file progression towards the apex.

When the last green bar is reached, a solid tone is emitted. The indication of the last green bar on the **X-Smart® Apex Locator** screen relates to the minor apical foramen file position – see Fig. 8.

### 6.6.3 Over-Instrumentation

A red “over” segment and an audio warning signal (rapid intermittent signal) indicate that the file has passed the apex – see Fig. 8.

### 6.6.4 Completion of the measurements

- Before unplugging the measurement cable from the device receptacle, disconnect the lip clip and the file clip from the patient.
- Move the file stopper to the selected reference point on the tooth.
- Gently remove the file from the canal and measure the apical length between the stopper and the file tip.

---

#### Note

Determination of working length for canal shaping is a subject of dentist’s professional judgment. In most cases subtraction of 0.5 mm from the measured apical length provides clinically acceptable working length. Nevertheless, in each case the dentist should define proper working length based on his experience, apex locator readings, radiographs and other available data.

---

## 6.7 Sound Adjustment

**X-Smart® Apex Locator** is equipped with an audio signal which enables monitoring of the progression of the file within the canal in addition to visual monitoring.

The volume can be adjusted by pressing the  Volume icon on the bottom left corner of the screen and sliding the scroller.


Fig. 9

## 6.8 Device Settings

To enter the Settings Menu, tap the  Settings icon on the top left corner on the screen. Tap on the required feature to activate or adjust its setting. Use the Up  and Down  arrows icons to see all menu features.


Fig. 10


Fig. 11

Setting Menu Features:

Dr's Choice	Set optional Dr's Choice apical line in the Apical region
Language	Set the user interface language
Check Mode	Check the operation of both apex locator/cables
Demo	Activate for demo purpose
Brightness	Adjust the display brightness
Sleep Mode	Set the operation time of the display before switching it automatically off
Auto Shutdown	Set the operation time of the device before switching it automatically off
 ECO Mode	Automated battery saving mode
About	Device firmware characteristics

## 6.8.1 Dr's Choice

This feature enables to mark an individual predetermined reference position at the required distance from the apex. This variable apical line can be set between the first green bar and the last green bar. When Dr's Choice apical line is set, clear visual and audio indication is given that the file tip has reached this pre-selected position. To set the Dr's Choice apical line or to modify the apical line position, follow the next steps:

- Enter the Settings Menu and tap the Dr's Choice tab – see Fig 12. "Dr's Choice" will turn orange while activated and the Dr's Choice apical line will appear on the screen – see Fig. 13.


Fig. 12


Fig. 13

- Tap on the Dr's Choice apical line  and set its position anywhere within the green bars. After setting its position tap somewhere on the screen.
- To deactivate Dr's Choice apical line, go to settings menu and turn Dr's Choice off by tapping the Dr's Choice tab – see Fig. 12.

### Note

When the apex position is reached (last green bar), a solid tone is sounded, even if the Dr's Choice feature is activated.

## 6.8.2 Language selection

To modify the user interface language:

- Enter the Settings Menu and tap the Language tab.
- Select your preferred language from the available languages list – see Fig. 14.
- Use the Up ▲ and Down ▼ arrows icons to see all the available languages.


Fig. 14

## 6.8.3 Check Mode

In case measurement indication is not detected this built-in Check feature enables to automatically test the basic functions of the device in a first step, followed by the test of its accessories in a second step.

To activate the Check mode feature, follow the next steps:

- Disconnect the measurement cable/charger from the device.
- Enter the Settings Menu and select the Check Mode feature.
- An automated internal test will start, and the results will be shown on the screen. Either ok -  to indicate that the device fully operates or “Error – Please contact your distributor”.

### Note

If the device functional check is ok, proceed to the cables functional check.

- Connect the file clip and the lip clip to the measurement cable and plug the measurement cable to the device.
- Tap on the “Run File Clip check” button and follow the instructions on the screen.
- The results will be shown on the screen. Either ok -  or “Error” will be shown on the screen.

---

### Warning

Error message indicates that the accessories are not functioning properly (e.g. cable breakage), or the contact area is dirty. Please check the connections or change the file clip and try again.

For assistance, please contact your distributor.

---

#### 6.8.4 Demo Mode

The built-in demo mode is available to become acquainted with the device and to demonstrate its operation.

- Disconnect the charger from the device, if connected, and turn the device on.
- Enter the Settings Menu and select the Demo Mode.
- During the demo cycle the device operating sequence is shown on the screen. It can be paused and resumed by tapping the pause/play button on the screen. The demo cycles are repeated automatically until interrupted by the operator.
- To exit the demo mode, press the back  icon located on the top left corner of the screen.

#### 6.8.5 Sleep Mode


The Sleep Mode is a power saving state. The **X-Smart® Apex Locator** will switch to a Sleep Mode when the device is in non-use state. To adjust the time duration before entering the Sleep Mode, adjust the slider position between 1 min to 5 min.

To wake-up the **X-Smart® Apex Locator** from the Sleep Mode, simply tap anywhere on the display, connect lip clip to file clip or start your measurements.

### 6.8.6 Automatic Shutdown

Automatic shutdown of the **X-Smart® Apex Locator** can be adjusted between 5 min to 15 min of device non-use. Simply adjust the Automatic Shutdown slider to your preferred time duration.

To prolong battery life, it is recommended to switch the device off after usage by pressing the “ON / OFF” button.

### 6.8.7 ECO Mode

The ECO mode will automatically manage the power consumption of the **X-Smart® Apex Locator** to prolong the battery usage. When ECO mode is activated the Brightness, Sleep Mode and Automatic Shutdown features settings will be disabled.

# 7 Maintenance, Cleaning, Disinfection and Sterilization

## 7.1 General Recommendations

- The device does not contain user serviceable parts. The service and repair should be provided by factory trained service personnel only.
- After each use, all objects that were in contact with infectious agents should be cleaned using single-use wipes or a soft cloth impregnated with a disinfecting and detergent solution (a bactericidal, fungicidal and aldehyde free solution), according to manufacturer's directions. We recommend using only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval). Use of chemical agents may cause damage to the equipment. To clean the device housing and measurement cable, we recommend single-use wipes (e.g. DENTIRO® Wipes). After cleaning, all disinfectant residues should be removed. No visible impurities or liquid should remain on the device housing after disinfection.
- In addition, the lip clip, the file clip and the touch probe must be sterilized between treatments. Please note that the device housing, the charger, and the measurement cable cannot be sterilized.
- Follow the “Disinfection and sterilization procedure” described in section 7.2 for the lip clip, the file clip and the touch probe.
- The user is responsible for the sterility of the lip clip, the file clip and the touch probe for the first cycle and each further usage.
- All damaged accessories should be discarded, while dirty accessories should be cleaned and sterilized per the procedure described in section 7.2.
- The maximum number of sterilization cycles is:
  - File clip: 200.
  - Lip clip: 200.
  - Touch probe: 200.

## 7.2 Disinfection and Sterilization Procedure for lip clip, file clip and touch probe

### Foreword

For hygiene and sanitary safety purposes, the lip clip, the file clip and the touch probe must be cleaned, disinfected and sterilized before each usage to prevent any cross-contamination between patients. This concerns the first use as well as the subsequent uses.

### General recommendations

- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents.
- For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.
- Limitations and restrictions on reprocessing: the appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of color, are indications that the devices are not able to fulfil the intended use with the required safety level.
- Use only clean water in all cleaning and rinsing steps.

### Step-by-step procedure

	Operation	Instructions	Details and Warning
1	Preparation at the point of use prior to processing	<ul style="list-style-type: none"> <li>• No particular requirements</li> </ul>	
2	Preparation for decontamination/ preparation before cleaning	<ul style="list-style-type: none"> <li>• No particular requirements</li> </ul>	

	<b>Operation</b>	<b>Instructions</b>	<b>Details and Warning</b>
3	Cleaning: Automated	<ul style="list-style-type: none"> <li>The accessories are not intended for automated cleaning</li> </ul>	
4	Cleaning: Manual	<ul style="list-style-type: none"> <li>Clean the accessories with an adequate brush or towel soaked in a disinfectant solution</li> </ul>	<ul style="list-style-type: none"> <li>The file clip should be activated during cleaning process (pressed and released several times)</li> <li>After cleaning no visible impurities should remain on the accessories</li> </ul>
5	Disinfection	<ul style="list-style-type: none"> <li>Soak the required accessories in a disinfectant solution combined with proteolytic enzyme if possible.</li> <li>Rinse well the accessories in flowing water</li> </ul>	<ul style="list-style-type: none"> <li>Follow instructions given by the disinfectant manufacturer (concentration, immersion time, etc.)</li> </ul>
6	Drying	<ul style="list-style-type: none"> <li>No particular requirements</li> </ul>	
7	Maintenance, inspection and testing of the accessories	<ul style="list-style-type: none"> <li>No particular requirements</li> </ul>	
8	Packaging	<ul style="list-style-type: none"> <li>Pack the accessories in sterilization pouches</li> </ul>	<ul style="list-style-type: none"> <li>Check the validity period of the pouch given by the manufacturer to determine the shelf life of the sterilized items</li> <li>Use packaging which is resistant up to a temperature of 141°C (286°F)</li> </ul>

	<b>Operation</b>	<b>Operating mode</b>	<b>Warning</b>
9	Sterilization	<ul style="list-style-type: none"> <li>• Steam sterilization at 135°C (275°F) during 10 minutes in gravity type autoclave (Tabletop, N type)</li> <li>• Drying time after sterilization – 30 minutes</li> </ul>	<ul style="list-style-type: none"> <li>• Follow maintenance and operation procedures of the autoclave provided by the manufacturer</li> <li>• The only sterilization parameters to be used are those that have been validated and provided to the user in this User Manual</li> </ul>
10	Storage	<ul style="list-style-type: none"> <li>• Keep sterilized accessories in sterilization pouches in a dry and clean environment</li> </ul>	<ul style="list-style-type: none"> <li>• Sterility cannot be guaranteed if packaging is open or damaged (check the packaging before using the accessories)</li> </ul>

## 8 Troubleshooting

Please review the checklist below should you experience a problem with your **X-Smart® Apex Locator**. If the problem persists after following the proposed solutions, please contact your distributor.

### Warning

The following patient related factors may prevent accurate readings:

- Blocked root canals;
- Teeth with large apices;
- Root fracture or perforation;
- Metal crowns or bridges if they come into contact with the file or the lip clip.

#	Problem	Possible cause	Solution
1	The device does not turn on by pressing the "ON / OFF" button.	The button is functioning incorrectly.	Try pressing the "ON / OFF" button several times. If the device still does not turn on, please contact your distributor.
		The battery is discharged.	Charge the battery.
		Electronic malfunction.	Contact your distributor.
2	The device shuts off during the procedure.	The battery is low.	Charge the battery.
3	No sound during the procedure.	The sound control is set at "Mute" level.	Adjust the sound level by pressing the volume  icon.
4	Display is not steady during the procedure.	There is not a good contact between the lip clip and the oral mucosa.	Ensure a good contact between the mucosa and lip clip (Place the lip clip in the labial angle opposite the tooth to be treated).
		The file clip is soiled.	Clean the file clip with a single-use wipes (e.g. DENTIRO Wipes)
		Deep caries provides a conductive path outside the canal.	Block the external conductive path.

#	Problem	Possible cause	Solution
		Perforation.	Remove the file, close the perforation and repeat the apex detection procedure, carefully inserting the file into canal.
		Large lateral canal.	Try continuing the procedure by gently advancing the file.
5	The transmission of the electric signal is interrupted. The device does not show file progression inside the canal.	Bad electrical contact.	Perform the Check Mode sequence as described in the user manual, section 6.8.3
		The file clip is not properly connected to the file.	Place the file clip on the metal part of the file below the plastic handle.
		The root canal is obliterated.	Check the comparative X-ray image for hints.
		In the case of retreatment: old filling material residues may block the root canal.	Remove old root filling material residues prior to use.
		The root canal may be blocked by the remnants of a medication (e.g. calcium hydroxide).	Completely remove the remnants prior to use.
		Root canal is extremely dry.	Rinse the root canal with NaCl solution. Dry the access cavity with a cotton pellet/ air-blower.
		The selected file is too small for a large root canal.	If there is no parietal contact use larger ISO size file. Important: exactly fitting files lead to precise results.

#	Problem	Possible cause	Solution
		Electronic malfunction.	Contact your distributor.
6	Display reaction is erratic: apical foramen position or “over” appear on the screen before the apical region is reached.	Short circuit due to excess liquid (irrigation solution, saliva, blood) in the pulp chamber.	Dry the access cavity with a cotton pellet / air-blower. In case of excess bleeding wait until it has stops.
		A direct contact of the file with the gingiva or gingival proliferations, e.g. a fractured metal crown.	For isolation: <ul style="list-style-type: none"> <li>adequate preparation filling of access cavity;</li> <li>use a rubber dam.</li> </ul>
		A direct contact of the file with metal restorations (crown, parapulpal post, amalgam filling).	Isolate the file by inserting it into a small polyvinyl tube before use.

## 9 Warranty

**X-Smart® Apex Locator** is warranted for 24 months from the date of purchase. The accessories (cables, battery, etc.) are warranted for 6 months from the date of purchase. Within the warranty period the manufacturer undertakes, at its sole discretion, to repair or replace the faulty item without charge.

This product has been developed specifically for use in dentistry and is intended to be operated only by qualified dental professionals in accordance with the instructions contained in this manual. However, notwithstanding anything contained herein, the user shall at all times be solely responsible for determining the suitability of the product for the intended purpose and the method of its use. Any guidance on technology application offered by or on behalf of the manufacturer, whether written, verbal or by demonstration, shall not relieve the dental professional from his/her obligation to control the product and to make all professional judgments regarding its use.

Except for the warranties specifically set forth in this manual, the manufacturer provides no warranties or guarantees of any kind covering the product, expressed or implied, including, without limitation, any warranties as to merchantability or fitness for a particular purpose. Any claim for damage or breakage to the product in transit should be made to the carrier promptly upon discovery.

The warranty is valid for normal usage conditions. Any damage caused by accident, abuse, misuse, or as a result of service or modification other than by a person authorized by the manufacturer will render the warranty void.

## 10 Disclaimer

The manufacturer, its representatives and its distributors shall have no liability or responsibility to customers or any other person or entity with respect to any liability, loss or damage caused or alleged to be caused directly or indirectly by equipment sold or furnished by us, including, but not limited to, any interruption of service, loss of business or anticipatory profits, or consequential damages resulting from the use or operation of the equipment.

The manufacturer reserves the right to implement changes and modifications of the product at any time, to revise this publication and to make changes in the contents hereof without obligation to notify any person of such changes, modifications, or revisions.

## 11 Certification

**X-Smart® Apex Locator** complies with the following standards: IEC 60601-1 (Safety) and IEC 60601-1-2 (Electromagnetic compatibility), including conducted and radiated immunity tests as specified for equipment of Group 1 Class B.

**X-Smart® Apex Locator** is covered by the “CE Marking of conformity” certificate. The device bears the following CE identification mark:

## 12 European Authorized Representative

European Authorized Representative who has been empowered to enter into commitments in our behalf:

CEpartner4U B.V.

ESDOORNLAAN 13, 3951 DB MAARN, THE NETHERLANDS.

Contact Person: Mr. Ton Pennings

Tel. +31 343 442 524, Fax +31 343 442 162

Email: [office@cepartner4u.com](mailto:office@cepartner4u.com)

## 13 Disposal of the Product

Recycling: PLEASE DO NOT THROW AWAY! This product and all its components must be recycled through your supplier.

## 14 Reporting of Incident to Manufacturer

Users shall report to the Manufacturer any serious incident related to the device, using this email: [info@forumtec.net](mailto:info@forumtec.net). Any personal information contained in this email is considered confidential and is accessible only by authorized personnel.

Any serious incident that has occurred in relation to the device should be also reported to the competent authority of the Member State in which the user and/or patient is established.

## 15 Technical Characteristics

**X-Smart® Apex Locator** electronic apex locator is a programmable electrical medical device, belongs to the following category of medical devices:

- Internally powered equipment
- Type BF Applied Parts
- Not suitable for use in the presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide
- Continuous operation
- Expected service life: 3 years
- Ingress of liquids – not protected

- The device is intended for indoor use only
- Environmental conditions during storage/transportation:
  - Temperature: -20°C to +60°C (-4°F to 140°F)
  - Relative humidity: 10% to 90%, non-condensing
  - Atmospheric pressure: 106 kPa to 50 kPa.
- Environmental conditions during device usage:
  - Temperature: 10°C to +40°C (50°F to 104°F)
  - Relative humidity: 10% to 90%, non-condensing
  - Atmospheric pressure: 106 kPa to 70 kPa.

**X-Smart® Apex Locator** is intended for use in electromagnetic environment specified for equipment of Group 1 Class B.

**Specifications:**

Dimensions (W x D x H)	95 x 111 x 146 mm
Weight	500 gr.
Type of screen	Color TFT capacitive touch display
Supply	3.6V Lithium-ion rechargeable battery 3,200 mAh
Switching charger	Input: 100-240 V AC ~ 50-60 Hz Output: 5V DC, 2,000 mA

## 16 Identification of Symbols

Symbols used in these directions for use, packaging, device and parts.

Symbol	Identification
	Serial number
	Catalogue number
	Lot number
	Direct current (connection for power supply)
	Manufacturer
	Date of manufacture
	Class II equipment
	Type BF applied part
	Consult instructions for use
	Refer to instruction manual/booklet
	Recycling PLEASE DO NOT THROW AWAY! This product and all its components must absolutely be recycled through your distributor
	Temperature limitation

	Humidity limitation
	Atmospheric pressure limitation
	Additional information, explanation on operation and performance
	Warning
	Medical Device
	Authorized Representative in the European Community
	Authorized Representative in Switzerland
	Authorized Representative in the United Kingdom
	UK Conformity Assessed marking
	INMETRO marking
	CE marking
	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner
	Sterilizable in a steam sterilizer (autoclave) at temperature specified
	Non sterile
	Indicates the number of pieces in the package. "X" is replaced by the number of pieces in the package
	Keep Dry

# Appendix

## Electromagnetic Compatibility

### Notes:

- **X-Smart® Apex Locator** requires special precautions with regard to electromagnetic compatibility.
- It must be installed and prepared for use as described in section 6.6.1 “Getting Started”.
- Certain types of RF wireless communication equipment such as mobile telephones are likely interfere to with **X-Smart® Apex Locator**.
- The recommended radiation levels of RF wireless communication equipment specified in this paragraph must therefore be complied with.
- **X-Smart® Apex Locator** must not be used near or on top of another device. If this cannot be avoided, it is necessary – before clinical use – to check the equipment for correct operation under the conditions of use.

## Electromagnetic Emissions

### Notes:

- **X-Smart® Apex Locator** is intended for use in the professional healthcare facility electromagnetic environment specified in the tables below.
- The user and/or installer of the unit must ensure that it is used in such an environment.
- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

<b>Declaration - electromagnetic emissions</b>		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1 Class A	The <b>X-Smart® Apex Locator</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	Class A	The <b>X-Smart® Apex Locator</b> is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the <b>X-Smart® Apex Locator</b> or shielding the location.
Voltage Fluctuations and Flicker IEC 61000-3-3:2013	Complies	

<b>Declaration – Electromagnetic Immunity</b>			
<b>Immunity test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	8 kV contact 2, 4, 8, 15 kV air	8 kV contact 2, 4, 8, 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transients / bursts, IEC 61000-4-4	2 kV for power supply lines 1 kV for input/ output lines	2 kV for power supply lines Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5	1 kV Line(s) to line(s) 2 kV Line(s) to earth 2kV Signal (input/output) to earth	1 kV Line(s) to line(s) 2 kV Line(s) to earth N/A	Mains power quality should be that of a typical commercial or hospital environment.

**Declaration – Electromagnetic Immunity**

<b>Immunity test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
<p>Voltage dips, short interruptions and voltage variations on power supply input lines</p> <p>IEC 61000-4-11</p>	<p>0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</p> <p>0% UT; 1 cycle and 70% UT; 25/30 cycles</p> <p>Single phase at 0° 0% UT; 250/300 cycles</p>	<p>0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</p> <p>0% UT; 1 cycle and 70% UT; 25/30 cycles</p> <p>Single phase at 0° 0% UT; 250/300 cycles</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the <b>X-Smart® Apex Locator</b> requires continued operation during power mains interruptions, it is recommended that the <b>X-Smart® Apex Locator</b> be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50/60 Hz) magnetic field</p> <p>IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

**Declaration – Electromagnetic Immunity**

IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V, 6 V</p> <p>3 V/m</p> <p>3 V/m from 0.15 to 80 MHz; 6 V/m from 0.15 to 80 MHz and 80% AM at 1 kHz</p> <p>3 V/m from 80 MHz to 2.7 GHz</p>	<p>3 Vrms, 6 V</p> <p>3 V/m</p> <p>3 V/m from 0.15 to 80 MHz; 6 V/m from 0.15 to 80 MHz and 80% AM at 1 kHz</p> <p>3 V/m from 80 MHz to 2.7 GHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the <b>X-Smart® Apex Locator</b>, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$ $d = \left[\frac{12}{V_2}\right]\sqrt{P}$ $d = \left[\frac{12}{E_1}\right]\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{23}{E_1}\right]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

**Recommended separation distances between portable and mobile RF communications equipment and the X-Smart® Apex Locator**

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz outside ISM bands	150 kHz to 80 MHz in ISM bands	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$	$d = \left[\frac{12}{V_2}\right]\sqrt{P}$	$d = \left[\frac{12}{E_1}\right]\sqrt{P}$	$d = \left[\frac{23}{E_1}\right]\sqrt{P}$
0.01	0.12	0.2	0.4	1
0.1	0.37	0.64	1.3	2.6
1	1.17	2	4	8
10	3.7	6.4	13	26
100	11.7	20	40	80

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Maximum power (W)	Distance (m)	Immunity test level (V/m)	Compliance level (V/m)
385	380 – 390	TETRA 400	Pulse modulation <sup>b)</sup> 18 Hz	1.8	0.3	27	27
450	430 – 470	GMRS 460, FRS 460	FM <sup>c)</sup> ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
710	704 – 787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217 Hz	0.2	0.3	9	9
745							
780							
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18 Hz	2	0.3	28	28
870							
930							
1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation <sup>b)</sup> 217 Hz	2	0.3	28	28
1845							
1970							
2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217 Hz	2	0.3	28	28
5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217 Hz	0.2	0.3	9	9
5500							
5785							

<sup>a)</sup> For some services, only the uplink frequencies are included.

<sup>b)</sup> The carrier shall be modulated using a 50 % duty cycle square wave signal.

<sup>c)</sup> As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

**PAGE INTENTIONALLY LEFT BLANK**

The latest revision can be accessed through this email:  
info@forumtec.net

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Technologies (96) Ltd.**